

SMART TRAC II

RAPID ANALYSIS SYSTEM

Analyze Moisture/Solids,
Fat and Oils in:

Dairy Products


Meat Products

Snack Foods

Condiments

Soaps

and more...




ACCURATE FAT AND MOISTURE ANALYSIS HAS NEVER BEEN EASIER!

Analyze for moisture/solids with SMART Turbo



Roll sample and place in Trac Tube



Insert into NMR chamber for analysis



Precise results in 8 seconds!





FAST, ACCURATE ANALYSIS for MOISTURE/SOLIDS, FATS and OILS

A Winning Combination!



RAPID ANALYSIS SYSTEM
For Moisture/Solids, Fats and Oils



The SMART Trac II partners advanced, high performance NMR technology with the fast, reliable drying capability of the SMART Turbo Microwave Moisture/Solids Analyzer to create a system which yields fast, accurate results without solvents or calibration. Results for the percentage of moisture, solids and fat are simple and easy.


- Analyze twice as many samples per day
- Green! No solvents, 95% less energy needed
- No costly calibration maintenance
- More accurate than FTIR and NIR
- Official AOAC Method 2008.06

SMART Turbo Microwave Moisture/Solids Analyzer

The SMART Turbo provides rapid, accurate moisture/solids analysis for a broad range of products. The system's patented IntelliTemp™ temperature control ensures that all moisture is completely removed without overheating the sample for better precision. SMART Turbo accurately analyzes a wide variety of sample types including solids, powders, slurries and liquids (0.1 - 99.9% moisture).

High Performance NMR Technology

NMR directly measures fats and oils. Unlike indirect methods, which only measure the surface of a sample, NMR detects fats and oils throughout the sample, yielding an accurate analysis that is independent of sample uniformity and not affected by changes in color or texture.



Easy Method Transfer

SMART Trac II has standard methods that are included with the instrument. Simply select from a list of methods and start analyzing without additional calibration. Custom methods are also available. Contact our applications support group and we can e-mail them to you directly. Combine this with the incredible stability of the SMART Trac II and you have a remarkably easy-to-use instrument that you never have to re-calibrate.

Once a method is developed on a PARENT; it can be copied to a CHILD and used immediately. You get the same results whether you are transferring methods across the lab or across the ocean.



RAPID ANALYSIS SYSTEM

For Moisture/Solids, Fats and Oils



SYSTEM STATUS AT A GLANCE

Ready 

Analyzing 

RAPID ANALYSIS
newly designed NMR
analyzes fat in 8 seconds

NEVER RECALIBRATE
all operators and all shifts get
the same results

COMPACT NMR
state of the art
neodymium magnet

SMALLER SIZE/LIGHTER WEIGHT
fits easily in laboratory
environments

WEIGHT STORED AUTOMATICALLY
by a built-in analytical balance that
continuously weighs the sample
during the drying process



GREEN METHOD
No solvents and SMART Trac II uses 95% less energy than
conventional methods

RAPID ANALYSIS
Turbo Technology for the fastest
moisture/solids analysis

EXCELLENT PRECISION
with proprietary Intelli-Temp™
temperature feedback system

ON-BOARD PRINTER

LIMS READY
ethernet and serial ports

TEST PRECISION TO 0.01 %
Results in % moisture, % solids & % fat
Also calculates %protein


STORES 100 METHODS & 300 TEST RESULTS

SIMPLE USER INTERFACE

PROGRAMMING IN 5 LANGUAGES
(English, French, Spanish, Italian & German)

EASY ACCESS
to microwave chamber and balance pan

*Only CEM Has The
IntelliTemp™ Advantage*




A BETTER SOLUTION

WHAT IS NMR & HOW DOES IT WORK?

NMR quantifies fats and oils in samples and eliminates the need for time consuming extractions and solvents. Nuclear Magnetic Resonance (NMR) is the same technique as Magnetic Resonance Imaging (MRI), which has been widely used in the medical profession for years to image the human body.

CEM combines microwave drying with NMR to accurately measure fats and oils. The SMART Turbo removes all of the water molecules from the sample, eliminating signal interference and allowing the NMR to accurately measure the amount of fat content in a wide variety of food products.



The word "nuclear" in the name of the technique refers to the instrument's ability to analyze the nuclei of the sample's atoms. NMR technology does not generate or use ionizing radiation.

Rapid Moisture/Solids Fat Analysis

ANALYZE TWICE AS MANY SAMPLES EVERY DAY		
(min:sec)		
	SMART Trac I	SMART Trac II
Milk	4:45	2:51
Cream	3:36	2:24
Cheddar	3:59	2:03
Ketchup	3:52	2:23
Hot Dog	6:00	2:21
Beef	5:02	3:30



NEW SMALLER SIZE

In most laboratories, benchtop space is at a premium, so we designed the SMART Trac II with a new, compact magnet that is smaller than any comparable instrument. The system requires no external computer and the flat panel LCD display can be used freestanding or mounted on the back of the magnet. The magnet and electronics are housed in separate units so that the electronics unit can be placed on the floor, under the bench, or even on a shelf to save bench space.



DELIVERS THE SAME QUALITY PERFORMANCE UP TO 10 FEET AWAY

SYSTEM SPECIFICATIONS

MOISTURE/SOLIDS RANGE	0.01% to 99.99% in liquids, solids and slurries, 0.01% resolution
FAT RANGE	0.01% to 99.99% in liquids, solids and slurries, 0.01% resolution
BALANCE CAPACITY	100 grams, 0.1 mg readability
PROGRAM/DATA STORAGE	100 methods and 300 results
STANDARD SOFTWARE	Software with proprietary Turbo Technology and IntelliTemp™ Temperature Controlled Drying.
DATA ENTRY	Keypad with menu-driven software
DISPLAY	Black and white VGA (320 x 240) Flat panel LCD display
ACCESSORY PORTS	2 serial, RS 232, 9 pin ports for network connection, external balance, or bar code reader, Parallel Port for external printer, 4 USB ports, 1 Ethernet CAT 5 port
STANDARD PRINTER	Internal impact printer
NETWORK INTERFACES	RS232, USB, Ethernet
INSTRUMENT DIMENSIONS:	
SMART Turbo	22 in (w) x 23 in (d) x 15 in (h) 56 cm (w) x 59 cm (d) x 37 cm (h)
SMART Trac II Magnet	14 in (w) x 14 in (d) x 13 in (h) 36 cm (w) x 36 cm (d) x 33 cm (h)
SMART Trac II Processor	12 in (w) x 13 in (d) x 16 in (h) 30 cm (w) x 33 cm (d) x 41 cm (h) May be located up to 9'8" (3 m) away from the magnet.
WEIGHT:	
SMART Turbo	55 lbs, 25 kg
SMART Trac II Magnet	110 lbs, 50 kg
SMART Trac II Processor	37 lbs, 17 kg
AMBIENT AIR TEMPERATURE	Temperature must be between 15 °C and 30 °C. For optimum stability and performance, the ambient temperature should not vary more than 5 °C/day.
VOLTAGE	100-240 V (50-60 Hz)
POWER	800 VA maximum, UPS recommended if power is not reliable
CURRENT	
SMART Turbo	10 amps (110-127 V, 60 Hz); 5 amps (220-240 V, 50/60 Hz) 10 amps (100 V, 50/60 Hz); 5 amps (200-208 V, 50/60 Hz)
SMART Trac II Magnet & Processor Power Conditioner Requirements	6.3 amps (100-240 V, 50/60 Hz) 2.1 kVA 120 V, 60 Hz 1.8 kVA 240 V, 50 Hz
RF PULSE GENERATOR:	Pulse power 300 W nominal. Pulse times variable in 100 ns increments. Transmit and receive phases selectable 0, 90, 180 and 270°. Nominal 90° pulse times 3 µs (10 mm probe) and 4.5 µs (18 mm probe)
MAGNET:	Permanent, thermally stabilized, 0.55 T (23.4 MHz proton), Homogeneity better than 10 ppm
SIGNAL DETECTION:	Dual channel (quadrature) detection with programmable low-pass filtering. Programmable data acquisition rate up to 10 MHz per pair of points.
PATENTS	6,548,303 6,548,304 7,125,721 7,220,591



Winner of the R&D 100 and the NSPE New Product Awards



CEM has been an ISO-certified facility since 1994.

CEM Corporation
PO Box 200
Matthews, NC 28106
Tel: 800-726-3331
Tel: 704-821-7015
Fax: 704-821-7894
Email: info@cem.com
www.cem.com

France
 CEM µWave S.A.S.
 Immeuble Ariane
 Domaine Technologique de Saclay
 4, rue Rene' Razel
 91892 ORSAY Cedex
 Tel: 33 (01) 69 35 57 80
 Fax: 33 (01) 60 19 64 91
 Email: info.fr@cem.com
 www.cemfrance.fr

Germany, Austria, Switzerland
 CEM GmbH
 Carl-Friedrich-Gauss-Str.9
 47475 Kamp-Lintfort
 Tel: (49) 2842-9644-0
 Fax: (49) 2842-9644-11
 Email: info@cem.de
 www.cem.de

Ireland
 CEM Technology (Ireland) Ltd.
 Sky Business Centre
 9a Plato Business Park
 Damastown
 Dublin 15
 Tel: +353 (0) 1 885 1752
 Fax: +353 (0) 1 885 1601
 Email: info.ireland@cem.com
 www.cemmicrowave.co.uk

Italy
 CEM S.R.L.
 Via Dell' Artigianato, 6/8
 24055 Cologno al Serio (Bg)
 Tel: (39) 35-896224
 Fax: (39) 35-891661
 Email: info.srl@cem.com
 www.cemmicroonde.com

Japan
 CEM Japan K.K.
 2-18-10 Takanawa
 Minato-ku, Tokyo
 108-0074
 Tel: +81-3-5793-8542
 Fax: +81-3-5793-8543
 Email: info@cemjapan.jp
 www.cemjapan.co.jp

United Kingdom
 CEM Microwave Technology Ltd.
 2 Middle Slade
 Buckingham Industrial Estate
 Buckingham MK181WA
 Tel: (44) 1280-822873
 Fax: (44) 1280-822873
 Email: info.uk@cem.com
 www.cemmicrowave.co.uk